

**BASIC**  
SHOWER TROLLEY

ARJO



# BASIC SHOWER TROLLEY

The ARJO *Basic Shower Trolley* is an excellent aid when showering recumbent residents. The *Basic Shower Trolley* allows you to transfer your resident directly from the bed to the shower and back again, with no manual lifting.



Equipped with a head pillow and foldable side supports with safety lock, the ARJO *Basic Shower Trolley* is height adjustable, helping to reduce postural load.



The *Shower Panel*, available in three different sink alternatives, provides an efficiently controlled delivery system for assisted showering. The unit features an accurate thermostatic mixer with shut-off. As well as the resident shower, panels are fitted with an equipment cleaning shower to ensure good system hygiene routines.

# ARJO

...with people in mind.

ARJO, Inc.  
50 N. Gary Avenue  
Roselle, Illinois 60172  
Phone: 800-323-1245  
Fax: 888-594-2756  
[www.arjo.com](http://www.arjo.com)

The soft mattress is comfortable to lie on and the wide and sturdy construction of the trolley enables easy access when showering the resident. The trolley is easily adjusted from either side to a comfortable working height.



The ARJO *Maxi Slide* or *Maxi Tube*, when used in conjunction with the *Basic Shower Trolley*, simplifies the transfer from bed.

## Product Information

Height adjustable between 22.8" (570 mm) and 35.6" (890 mm)

All four castors have brakes

Equipped with head pillow and foldable side supports with safety lock

Length 76" (1900 mm)

Width 30.4" (760 mm)


Max load 330 lbs (150 kg)

## Accessory

Drain hose

## Resident Gallery

Suitable for Emma



**Emma**, who is almost completely bedridden and totally dependant

Please contact ARJO for further information on the Resident Gallery™.

Use only ARJO Designed Parts, i.e. parts specifically designed for the purpose, on equipment and products supplied by ARJO. Our Conditions of Sale specify that ARJO is not liable for injuries attributable to the use of inadequate parts. Continuous development is part of our policy, and we therefore reserve the right to change specifications without notice.

Basic is a trademark belonging to ARJO International AG.



FM 24321  
ISO 9001 • BS 5750 Part 1 • EN 29001  
ARJO Hospital Equipment AB, Sweden

MEMBER OF THE GETINGE GROUP