

Bedside Cabinet

GETINGE GROUP

Durable and modern bedside furniture



The **ArjoHuntleigh Bedside Cabinet** is a two-drawer bedside cabinet offering storage space for the patient's belongings.


- Wipe-clean top surface and drawers
- Towel/bag hook
- Lockable top drawer
- Large lower drawer with bottle holder
- Castors with kick-stop brakes
- Optional perimeter rail




Colours and options



Bedside cabinets are available in 4 colours and 3 wood finishes:



Specifications / Configurations	
Overall dimensions	Width 48.5 cm (19 inches), Depth 47.5 cm (18.5 inches), Height 74 cm (29 inches)
Product weight (approx.)	26 kg (57.3 lb)
Castor size	60 mm (6 cm / 2.3 inches)
Top surface	48 x 48 cm (18.9 x 18.9 inches)
Top drawer internal dimensions	Width 37.2 cm (14.5 inches), Depth 37.3 cm (14.6 inches), Height 8.6 cm (3.1 inches)
Top drawer front panel dimensions	Width 398 mm (15.7 inches), Thickness 8 mm (0.3 inches), Height 167 mm (6.6 inches)
Top drawer volume	20810 cm ³ (1270 in ³)
Bottom drawer internal dimensions	Width 37.2 cm (14.5 inches), Depth 37.3 cm (14.6 inches), Height 8.6 cm (3.1 inches)
Bottom drawer front panel dimensions	Width 398 mm (15.7 inches), Thickness 8 mm (0.3 inches), Height 348 mm (13.7 inches)
Bottom drawer volume	47175 cm ³ (2879 in ³)
Safe working load	 Bedside Cabinet: 100 kg (220 lb) evenly distributed Single Drawer: 10 kg (22 lb) evenly distributed
Flame rating	Classification M1 according to Standard NF P 92-507:2004
Material	
Panels	High Pressure Compact Laminate (HPCL)
Posts	Aluminium Extrusion Anodized Finish
Drawers	Injection Moulded Plastic
Material properties of HPCL	Resistance to scratches and impact from daily use Resistance to moisture, steam and heat

Available Configurations	
Colours	AB – AH Blue
	LL – Lilac
	LB – Light Blue
	LG – Light Grey
	WO – White Oak
	CB – Copper Beech
Options	NW – North Walnut
	BC0 Bedside Cabinet, no rails BC1 Bedside Cabinet, with rails

GETINGE GROUP

Getinge Group
 P.O Box 8861 SE-402 72
 Göteborg Sweden
www.getingegroup.com

MAQUET
 GETINGE GROUP

ARJOHUNTLEIGH
 GETINGE GROUP

GETINGE
 GETINGE GROUP

Getinge Group is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of ArjoHuntleigh, Getinge and Maquet. ArjoHuntleigh focuses on patient mobility and wound management solutions. Getinge provides solutions for infection control within healthcare and contamination prevention within life sciences. Maquet specializes in solutions, therapies and products for surgical interventions, interventional cardiology and intensive care.